

## Liner Sampler Sets



**Liner Sampler Sets** for collecting undisturbed samples for soil fertility, structure, volume weight, granular composition and technical determination in either soft or hard soils to a depth of 7m.

### Undisturbed samples collected in a liner for on site analysis

#### Advantages

- Liner allows easy removal of sample
- Sample can be viewed directly in the field
- Sample to a depth of 7m
- Can be hammered into position
- Sample can be analysed without touching
- Made from stainless steel
- Easily cleaned and decontaminated.

#### Disadvantages

- Not suitable for very stony soils.

#### Specifications

- Sample diameter 37mm (soft soils), 49mm (hard soils)
- Sample length 20cm (soft soils), 30cm (hard soils).

#### Alternatives

For top layers a Split Tube Sampler is a good alternative.